Good Quality Improves Life's Quality
PT. MERSIFARMA TM.
SUKABUMI - INDONESIA

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**TRANSFER PROTOCOL ANALYTICAL METHODS FOR
DETERMINING CONCENTRATIONS, RELATED
COMPOUND, AND PRESERVATIVE LEVELS BY
HPLC, EG AND DEG CONTAMINANTS BY GC-MS,
AND VISUAL IDENTIFICATION OF PIRACETAM IN
MERSITROPIL SYRUP**

Halaman	: 15 dari 24
Pages	: 15 of 24
No. Dokumen	: RD.02-PTP/102.00
Documen Number	:
Tanggal Dokumen Disetujui	:
Document Approved Date	:
Menggantikan No.	:
Replaces No.	:

- 7.8.6.2. Larutan Induk Uji kesesuaian sistem (Butilparaben 0,1 mg/mL; Propilparaben 0,1 mg/mL; Etilparaben 0,1 mg/mL; Metilparaben 0,1 mg/mL; dan asam p-hidroksibenzoat 0,1 mg/mL).

Main system Suitability Test Solution (Butylparaben 0.1 mg/mL; Propylparaben 0.1 mg/mL; Ethylparaben 0.1 mg/mL; Methylparaben 0.1 mg/mL; and p-hydroxybenzoic acid 0.1 mg/mL).

- Timbang seksama masing-masing 5 mg** Butilparaben, Propilparaben, Etilparaben, Metilparaben, dan asam p-hidroksibenzoat ke dalam labu tentukur masing-masing 50 mL.
Weigh carefully each 5 mg Butylparaben, Propylparaben, Etylparaben, Methylparaben, and p-hydroxybenzoate acid put each into 50 mL Volumetric Flask
- Tambahkan masing-masing 25 mL pelarut, sonikasikan 10 menit.
Add each 25 mL solvent, sonicate for 10 minutes.
- Tambahkan pelarut sampai tanda batas, kocok sampai homogen.
Add solvent up to volume, shake until homogenous.

- 7.8.6.3. Larutan Uji kesesuaian sistem (Butilparaben 0,01 mg/mL; Propilparaben 0,01 mg/mL; Etilparaben 0,01 mg/mL; Metilparaben 0,01 mg/mL; dan asam p-hidroksibenzoat 0,01 mg/mL).

System Suitability Test Solution (Butylparaben 0.01 mg/mL; Propylparaben 0.01 mg/mL; Ethylparaben 0.1 mg/mL; Methylparaben 0.01 mg/mL; and p-hydroxybenzoic acid 0.01 mg/mL).

System Suitability Test Solution (Consentration 0.01 mg/mL)

- Pipet masing-masing 1 mL Butilparaben, Propilparaben, Etilparaben, Metilparaben, dan asam p-hidroksibenzoat ke dalam labu tentukur 10 mL.
pipette each 1 mL Butylparaben, Propylparaben, Etylparaben, and p-hydroxybenzoate acid put into 10 mL Volumetric Flask.
- Tambahkan pelarut sampai tanda batas, kocok sampai homogen.
Add solvent, shake until homogenous.
- Saring larutan menggunakan syringe filter 0,45 µm PTFE.
Filter the solution using PTFE 0.45 µm syringe filter.
- Masukkan masing-masing konsentrasi ke dalam vial HPLC.
Put each concentration into HPLC vial.

- 7.8.6.4. Larutan Baku Internal (Etilparaben 0,013 mg/mL)

Internal Standard Solution (Ethylparaben 0.013 mg/mL)

- Pipet 2,6 mL larutan induk uji kesesuaian sistem Etilparaben ke dalam labu tentukur 20 mL.
Pipette 2.6 mL of Ethylparaben of main system suitability solution into a 20 mL volumetric flask.
- Tambahkan 10 mL pelarut, sonikasikan 10 menit.
Add 10 mL solvent, sonicate for 10 minutes.
- Tambahkan pelarut sampai tanda batas, kocok sampai homogen.
Add solvent, shake until homogenous.

Good Quality Improves Life's Quality

PT. MERSIFARMA TM.
SUKABUMI-INDONESIA

**PROTOKOL TRANSFER METODE ANALISA
PENETAPAN KADAR DAN CEMARAN ORGANIK
BAHAN BAKU PROPILPARABEN
TRANSFER PROTOCOL FOR ANALYSIS METHOD
FOR DETERMINING CONTENT AND ORGANIC
IMPURITY OF PROPYLPARABEN RAW
MATERIALS**

Halaman : 6 dari 11
Pages : 6 of 11
No. Dokumen : RD.02-PTB/090.00
Document Number :
Tanggal Dokumen Disetujui :
Document Approved Date :
Menggantikan No. :
Replaces No. :

FP std = Faktor pengenceran larutan standar (500)

Dilution factor of standard solution (500)

F Kadar Working Standar (WS) dalam persen (%)

Working Standard (WS) content in percent (%)

4.3. Penetapan cemaran organik

Determination of Organic Impurity

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Requirement In accordance with applicable product specifications.

4.3.2. Metode : High Performance Liquid Chromatography

Method High Performance Liquid Chromatography

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Run Time 2.5 times the retention time of Propylparaben

4.3.3.8. Waktu Retensi : 3.2 ± 2 menit

Retention Time 3.2 ± 2 minutes

4.3.3.9. End Time : 8.0 ± 2 menit

End Time 8.0 ± 2 minutes

4.3.3.10. Resolusi : ≥ 3.0 antara puncak etilparaben dan propilparaben (Larutan standar A)

≥ 3.0 between ethylparaben and propylparaben peaks (Standard solution A)

4.3.4. Larutan Standar A (Konsentrasi : 0.005 mg/mL)

Standard solution A (Concentration : 0.005 mg/mL)

4.3.4.1. Timbang seksama masing - masing 25 mg p-hydroxy benzoic acid, etilparaben RS dan Propylparaben WS.

Weigh accurately each 25 mg p-hydroxy benzoic acid, ethylparaben RS and Propylparaben WS.

4.3.4.2. Masukkan masing - masing ke dalam labu ukur 50 mL secara terpisah.


Put each into volumetric flask 50 mL separately

4.3.4.3. Tambahkan 20 mL pelarut, sonikasikan selama 10 menit.

Add 20 mL of solvent, sonicate for 10 minutes

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Add solvent up to volume, mix until homogeneous

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[↑](#) Bagikan

Ethylparaben

Pharmaceutical Secondary Standard; Certified Reference Material

★★★★★ (0) [Write a review](#) [Ask a question](#)**Sinonim:** Ethyl 4-hydroxybenzoate, EP, Ethylparaben[Masuk untuk Melihat Harga Organisasi & Kontrak.](#)

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