Good Quality Improve Life's Quality

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6/25/2026

9:52:21 AM

Item Number		RGN00055		Calcium carbonate QC		Buyer/Planner		Site		MTM	
Qty On Hand		0.0		UM Gr		Time Fence 0		Manufacturing Lead Time 0		MRP Required No	
Order Policy		POQ		Minimum Order 0		Safety Time 0		Pur/Mfg P Purchase LT 40		Mstr Sched Yes	
Order Period		0		Maximum Order 0		Safety Stk 0		Inspect No Inspect LT 0		Plan Orders Yes	
Order Qty		0		Ord Mult 0		Yield% 100.00%		Cumulative Lead Time 0			
Issue Policy		Direct						BOM/Formula Code			

**PROTOKOL TRANSFER METODE ANALISIS
KADAR KALSIUM, CEMARAN ORGANIK DAN
IDENTIFIKASI BAHAN BAKU KALSIUM PANTOTENAT
ANALYTICAL METHOD TRANSFER PROTOCOL
FOR ASSAY OF CALCIUM, ORGANIC IMPURITIES AND
IDENTIFICATION OF CALCIUM PANTOTHENATE RAW
MATERIALS**

Halaman	: 7 of 15
Pages	
No. Dokumen	: RD.02-PTB/127.00
Document No.	
Tanggal Dokumen Disetujui	:
Document Approval Date	
Menggantikan No.	:
Supersedes Number	

- 7.1.2.2. *Purified water*
Purified water
- 7.1.2.3. Natrium Hidroksida ex.Thermo (No.Cat 134070010)
Sodium Hydroxide ex.Thermo (No.Cat 134070010)
- 7.1.2.4. Dinatrium Etildiamintetraasetat ex.Kanto (No.Cat JIS K 8107)
Disodium Ethyldiaminetetraacetate ex.Kanto (No.Cat JIS K 8107)
- 7.1.2.5. Kalsium Karbonat ex.Sigma (No.Cat 239216)
Calcium Carbonate ex. Sigma (No.Cat 239216)
- 7.1.2.6. Biru Hidroksinaftol ex.Merck (No.Cat 1.04593.0025)
Hydroxynaphthol blue ex.Merck (No.Cat 1.04593.0025)
- 7.1.3. Pembuatan Larutan HCl 10%
Preparation of 10% HCl solution
- 7.1.3.1. Pipet 31.25 mL larutan HCl 32%, masukkan ke dalam labu ukur 100 mL.
Pipette 31.25 mL of 32% HCl solution, put into a 100 mL volumetric flask.
- 7.1.3.2. Larutkan dan encerkan dengan *purified water* hingga batas tanda.
Dissolve and dilute with purified water to the mark.
- 7.1.4. Pembuatan Larutan HCl 3 N
Preparation of 3 N HCl solution
- 7.1.4.1. Pipet 29.5 mL larutan HCl 32% ke dalam labu ukur 100 mL.
Pipette 29.5 mL of 32% HCl solution, put into a 100 mL volumetric flask.
- 7.1.4.2. Larutkan dan encerkan dengan *purified water* hingga batas tanda.
Dissolve and dilute with purified water to the mark.
- 7.1.5. Pembuatan Larutan NaOH 1 N
Preparation of 1 N NaOH solution
- 7.1.5.1. Timbang kurang lebih 8 gram NaOH, masukkan kedalam labu ukur 200 mL.
Weigh approximately 8 grams of NaOH, put it into a 200 mL measuring flask.
- 7.1.5.2. Larutkan dan encerkan dengan *purified water* hingga batas tanda.
Dissolve and dilute with purified water to the mark.
- 7.1.6. Pembuatan Larutan Dinatrium Edetat 0.05 M
Preparation of 0.05 M Disodium Edetat solution
- 7.1.6.1. Timbang kurang lebih 18.6 gram dinatrium etildiamintetraasetat, masukkan ke dalam labu ukur 1000 mL.
Weigh approximately 18.6 grams of disodium ethyldiaminetetraacetate, put it into a 1000 mL measuring flask.
- 7.1.6.2. Larutkan dan encerkan dengan *purified water* hingga batas tanda.
Dissolve and dilute with purified water to the mark.
- 7.1.7. **Pembakuan Larutan Titran Dinatrium Edetat 0.05 M**
Standardization of 0.05 M Disodium Edetat titrant solution
- 7.1.7.1. **Timbang seksama lebih kurang 200 mg kalsium karbonat** yang telah dikeringkan pada suhu 110°C selama 2 jam.
Accurately weigh approximately 200 mg of calcium carbonate, dried at 110°C for 2 hours.
- 7.1.7.2. Tambahkan 10 mL *purified water* hingga menjadi bubur.
Add 10 mL of purified water to form a slurry.
- 7.1.7.3. Tutup erlenmeyer.
Cap the erlenmeyer.
- 7.1.7.4. Tambahkan 2 mL larutan HCl 10%, kocok hingga larut dan homogen.
Add 2 mL of 10% HCl solution, shake until dissolved and homogeneous.
- 7.1.7.5. Tambahkan 100 mL *purified water*, kocok hingga homogen.
Add 100 mL of purified water, shake until homogeneous.
- 7.1.7.6. Titrasi dengan larutan dinatrium edetat 0.05 M sampai kurang lebih 30 mL.
Titrate with 0.05 M disodium edetat solution to approximately 30 mL.
- 7.1.7.7. Tambahkan 15 mL larutan NaOH 1 N dan 300 mg indikator biru hidroksinaftol.

**PROTOKOL TRANSFER METODE ANALISA
IDENTIFIKASI DAN PENETAPAN KADAR PADA
BAHAN BAKU **CALCIUM CARBONATE**
ANALYSIS METHOD TRANSFER PROTOCOL
FOR IDENTIFICATION AND ASSAY IN CALCIUM
CARBONATE RAW MATERIAL**

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Preparation of Sodium Hydroxide 1 N Solution

- Timbang seksama 4.0 g natrium hidroksida, masukkan ke labu ukur 100 mL.
Weigh carefully 4.0 g of sodium hydroxide, place in a 100 mL volumetric flask.
- Tambahkan 50 mL *purified water*, aduk sampai larut.
Add 50 mL of purified water and stir until dissolved.
- Tambahkan *purified water* sampai 100 mL, aduk sampai homogen.
Add purified water up to 100 mL, shake until homogeneous.

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Preparation of Dilute Hydrochloric Acid Solution

- Pipet 22.6 mL asam hidroklorida 32%, masukkan ke labu ukur 100 mL.
Pipette 22.6 mL of hydrochloric acid 32%, place in a 100 mL volumetric flask.
- Tambahkan 50 mL *purified water*, aduk sampai larut.
Add 50 mL of purified water and stir until dissolved.
- Tambahkan *purified water* sampai 100 mL, aduk sampai homogen.
Add purified water up to 100 mL, shake until homogeneous.

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Equivalence of EDTA with Calcium Carbonate

**PROTOKOL TRANSFER METODE ANALISA
IDENTIFIKASI, PENETAPAN KADAR, DAN
KAPASITAS PENETRALAN ASAM PADA
ANTASIDA DOEN TABLET
ANALYSIS METHOD TRANSFER PROTOCOL
FOR IDENTIFICATION, ASSAY, AND ACID
NEUTRALIZING CAPACITY IN ANTASIDA DOEN
TABLETS**

Halaman	: 12 dari 22
Pages	: 12 of 22
No. Dokumen	: RD.02-PTP/136.00
Documen Number	:
Tanggal Dokumen Disetujui	:
Document Approved Date	:
Menggantikan No.	: -
Replaces No.	:

7.5.11. Pembuatan larutan indikator *Eriochrome Black T*

Preparation of Eriochrome Black T indicator solution

7.5.11.1. Timbang seksama 200.0 mg *Eriochrome Black T*.
Weigh accurately 200.0 mg of Eriochrome Black T.

7.5.11.2. Tambahkan 15 mL *triethanolamine*.
Add 15 mL triethanolamine.

7.5.11.3. Tambahkan 5 mL *ethanol*, aduk hingga homogen.
Add 5 mL of ethanol, stir until homogeneous.

7.5.12. **Pembakuan larutan EDTA 0.05 M**

Standardization of 0.05 M EDTA solution

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7.5.12.2. *Weigh accurately 200 mg Calsium Carbonat (CaCO_3) which has been dried at 110°C for 2 hours, into a 250 mL Erlenmeyer flask.*

7.5.12.3. Tambahkan 10 mL *purified water*, goyangkan hingga menjadi bubur.
Add purified water 10 mL, shake it until it becomes porridge.

7.5.12.4. Tambahkan 2 mL HCl encer, kocok sampai homogen.
Add 2 mL HCl encer, shake until homogeneous.

7.5.12.5. Tambahkan 100 mL *purified water*, kocok sampai homogen.
Add purified water 100 mL, shake until homogeneous

7.5.12.6. Tambahkan 15 mL Natrium hidroksida LP, kocok sampai homogen.
Add 15 mL Sodium Hydroxide LP, shake until homogeneous.

7.5.12.7. Tambahkan 300 mg Indikator Hidroksinaftol P.
Add 300 mg of Hydroxynaphthol Indicator P.

7.5.12.8. Titrasi dengan larutan EDTA 0.05 M sampai titik akhir warna biru.
Titrate with 0.05 M EDTA solution until the blue end point.
Continue titration with 0.05 M EDTA solution until a blue end point occurs.

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Pipette 10.0 mL of 0.05 M EDTA solution, put it into a 250 mL Erlenmeyer flask.

7.5.13.2. Tambahkan 10.0 mL *Buffer acetic acid- ammonium acetate*.

Add 10.0 mL of acetic acid-ammonium acetate buffer.

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