Bagian : QC

Due Date : 15 - Sep - 26

Tanggal : 27 - Jul - 26

ANALISA PERMINTAAN BARANG

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No.	Desc.	No Cat.	Negara Asal	Qty (Pcs)	Price
1.	Camag Twin Trough Chamber for 20 x 20 cm plates, with stainless steel lid	022.5256	Switzerland	1	Rp. 14.255.000

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Mengetahui,

Lucky Pratama
Area Sales Manager

Good Quality Improves Life's Quality

PT. MERSIFARMA TM.

SUKABUMI-INDONESIA

**METODE ANALISA PRODUK
TRIFLUOPERAZINE HCl 5 TABLET
SALUT SELAPUT
PRODUCT ANALYSIS METHOD OF
TRIFLUOPERAZINE HCl 5 FILM
COATED TABLETS**

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Notes : Standard and samples solutions are protected from light

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Perform scanning on the sample solution and standard solution at wavelengths of 200-400 nm. The conditions for processing samples and standards are the same as for determining levels. Observe the wavelengths of the minimum and maximum absorbance obtained, with an acceptance criterion of ± 2 nm.

8.2. Uji B

Test B

- 8.2.1. Persyaratan : Sesuai dengan spesifikasi produk yang berlaku.
Requirement : In accordance with applicable product specifications.
- 8.2.2. Metode : Kromatografi Lapis Tipis
Method : Thin Layer Chromatography

8.2.3. Kondisi pengerjaan :
Working Conditions


- 8.2.3.1. Lempeng kromatografi : Silica Gel 60 F254
Plate chromatography : Silica Gel 60 F254
- 8.2.3.2. Fase Gerak : Aseton : Amonium Hidroksida (200 : 1)
Mobile Phase : Aceton : Ammonium Hydroxide (200 : 1)
- 8.2.3.3. Volume penotolan : 5 μ L
The volume of bottling : 5 μ L
- 8.2.3.4. Pelarut : Metanol
Solvent : Methanol

8.2.4. Pembuatan Larutan Standar (Trifluoperazin HCl 1,2 mg/mL)

Preparation of Standar Solution (Trifluoperazine HCl 1.0 mg/mL)

- 8.2.4.1. Timbang saksama 12 mg Trifluoperazin Hidroklorida WS, masukkan ke dalam labu tentukur 10 mL.

Weigh accurately 12 mg Trifluoperazine hydrochloride WS, put into a 10 mL volumetric flask.

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PT. MERSIFARMA TM
DIVISI MANUFACTURING

METODE ANALISA PRODUK
RANITIDINE 50 mg / 2 mL INJEKSI
PRODUCT ANALYSIS METHOD OF
RANITIDINE 50 mg / 2 mL INJECTION

DEPARTEMEN
RESEARCH AND DEVELOPMENT

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Calculate the recovery (%) of Ranitidine in Ranitidine 50 mg / 2 mL Injection according to the formula :

$$\% \text{ Ranitidine} = \frac{C}{B} \times A$$

Keterangan :

Description :

C = Bobot isi tiap ampul Ranitidine 50 mg / 2 mL injeksi (mg)

Weight the content of each Ranitidine 50 mg / 2 mL Injection Ampoule (mg)

B = Bobot rata-rata isi ampul Ranitidine 50 mg / 2 mL injeksi (mg)

Average weight of Ranitidine 50 mg / 2 mL Injection Ampoule (mg)

A = Hasil rata-rata penetapan kadar (%)

Average of assay result (%)

8.1.3. Perhitungan nilai penerimaan

Calculation of Acceptance Value

8.1.3.1. Hitung nilai keberterimaan dengan rumus:

Calculate the acceptance value using the following formula:

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Acceptance constant, depending on the number of units tested (for example, k = 2.4 if n = 10, k = 2.0 if n = 30)

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Requirement In accordance with the applicable finished drug specifications.

9.2. Metode : Kromatografi Lapis Tipis

Method Thin Layer Chromatography

PT. MERSIFARMA TM
SUKABUMI - INDONESIA

METODE ANALISA PRODUK
DIAZEPAM 5 TABLET
PRODUCT ANALYSIS METHOD
DIAZEPAM 5 TABLETS

DEPARTEMENT OF
RESEARCH AND DEVELOPMENT

SECTION OF
ANALYTICAL DEVELOPMENT

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Record the friability of each tablets.

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Requirement : *Meets the requirements for product specification*

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Put 1000 mL of purified water into the beaker glass of the disintegration tester.

7.3. Atur suhu air murni $37^{\circ}\text{C} \pm 2^{\circ}\text{C}$.

Set the pure water temperature to $37^{\circ}\text{C} \pm 2^{\circ}\text{C}$.

7.4. Masukkan Diazepam 5 tablet ke dalam masing-masing vessel.

Put Diazepam 5 tablets into each vessel.

7.5. Jalankan alat *disintegration tester*.

Put Diazepam 5 tablets into each vessel.

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Record the disintegration time of each tablet.

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Compare the relative retention time of the major peak of the sample solution with that of the standard solution obtained in the assay.

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

Acceptance criteria : In accordance with applicable product specifications.

5.2. Cara Pengujian :

Procedure :

5.2.1. Siapkan 10 Botol Sampel Paracetamol Syrup.


Prepare 10 sample bottles of Paracetamol Syrup.

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Transfer the contents of each bottle one by one into a 100 mL measuring cup.

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Note the volume stated on the measuring cup.

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Requirement : In accordance with applicable product specifications

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Compare the respons peak area of sample chromatogram to the standar solution in accordance with the applicable product specifications.

6.2. Identifikasi dengan metode KLT

Identification with Thin-Layer Chromatography method

6.2.1. Metode : Kromatografi Lapis Tipis

Method : Thin Layer Chromatography

6.2.2. Kondisi Pengerjaan KLT

TLC Operating Conditions

6.2.2.1. Lempeng kromatografi : KLT Silica Gel 60 F254

Absorber : F254 60 silica gel plate.

6.2.2.2. Pengamatan bercak : UV 254 nm

Spot observation : UV 254 nm

6.2.2.3. Jarak rambat eluen : 8 cm

Eluent distance : 8 cm

6.2.2.4. Eluen : Diklorometan : metanol (4:1)

Eluent : Dichloromethane : Methanol (4 : 1)

6.2.2.5. Volume penotolan : 10 µL

Spotting volume : 10 µL

6.2.2.6. Pelarut : Metanol

Solvent : Methanol

6.2.3. Larutan standar

Standard Solution

6.2.3.1. Timbang seksama 54 mg Parasetamol WS, masukkan ke dalam labu ukur 25 mL.

Weigh accurately 54.0 mg Paracetamol WS. put into 25 mL volumetric flask.

6.2.3.2. Tambahkan 15 mL pelarut, sonikasikan ± 10 menit.

Add 15 mL solvent, sonicate ± 10 minutes