ANALISA PERMINTAAN BARANG

No	Kode Barang	Nama Barang	Satuan	Stok Saat Ini (A)	Rencana Kebutuhan (B)	Kekurangan (A -B + C)	Pending Penerimaan (C)	Jumlah	Realisasi Pembelian	Realisasi Penggunaan
1	IBS00049	RS Fluphenazine Hydrochloride	mg	0	200	-200	0	200	0	0
2	IBS00469	RS Fluphenazine Decanoate Dihydrochloride	mg	0	100	-100	0	100	0	0

Catatan :

- 1. Fluphenazine Hidroklorida Ex. BPFI, No. Kontrol: B0120408 (1 Vial @200 mg)
- 2. Fluphenazine Decanoate Dihydrochloride Ex. USP Cat.no. 1280803 (1 Vial @100 mg)
- Digunakan untuk baku pembanding analisa produk Fluphenazine Decanoate Injeksi 25mg/mL.
- Pemesanan pertama kali Dept.QC

Latar belakang permintaan :

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Ready to ship

uspi

Fluphenazine Decanoate Dihydrochloride (100 mg)

Catalog No : 1280803

CAS RN® 2376-65-0

Molecular Formula: C32H44F3N3 O2S.2HCl

Product Type: Reference Standard

1

+

Add to cart

Add to Favorites

Shipping Information

Related Products

THIS ITEM HAS BEEN IDENTIFIED AS A POTENTIALLY DANGEROUS GOOD. SPECIAL SHIPPING METHODS AND CHARGES MAY APPLY.

\$329.00

Frequently Purchased

Tenivastatin Calcium (2 x 15 mg) (CO LD SHIPMENT REQUIRED) (Calcium (3R,5R)-7-((1S,2S,6R,8S,8aR)-8-((2,2-dimethylbutanoyl)oxy)-2,6-dimethyl-1,2,6,7,8,8a-hexahydronaphthalen-1-yl)-3,5-dihydroxyheptanoate (1:2) monohydrate)

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Showing 1 to 1 of 1 entries (filtered from 660 total entries)

Nomor	Nama Baku Uji	Keterangan	No. Kontrol	No. Kontrol Sebelumnya	Bobot	Ket. Pengujian
374	Flufenazin hidroklorida	BPFI	B0120408	-	200 mg	Kuantitatif

PreviousNext

Good Quality Improves Life's Quality

PT. MERSIFARMA TM.
SUKABUMI-INDONESIA

PROTOKOL TRANSFER METODE ANALISIS IDENTIFIKASI,
KADAR, BATAS FLUPHENAZINE DAN BATAS CEMARAN
YANG TERELUSI LAMBAT DALAM PRODUK JADI
FLUPHENAZINE DECANOATE INJEKSI 25 mg/mL
TRANSFER PROTOCOL OF ANALYSIS METHOD OF
IDENTIFICATION, ASSAY, LIMIT OF FLUPHENAZINE AND
LIMIT OF SLOWLY ELUTERATED IMPURITIES IN FINISHED
PRODUCT FLUPHENAZINE DECANOATE INJECTION 25
mg/mL

Halaman : 7 dari 15
Pages : 7 of 15
No. Dokumen : RD.02-PTP/005.00
Document No.
Tanggal Dokumen Disetujui :
Document Approved Date :
Menggantikan No. :
Supersedes Number :

No.	Nama Alat Tool Name	Jumlah Amount
11	Labu ukur 10 mL Volumetric flask 10 mL	6
12	Labu ukur 50 mL Volumetric flask 50 mL	2
13	Magnetik Stirer Magnetik Stirer	1

No.	Nama Alat Tool Name	Jumlah Amount
14	Saringan Vakum Vacuum Filter	1
15	Piknometer Piknometer	1
16	Pipet tetes Dropper pipet	1

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Materials Used

No.	Nama Bahan Material Name	Nomor Katalog Catalog Number	Manufature
1	Asetonitril Acetonitrile	A998-4	Thermo
2	Metanol Methanol	A452-4	Thermo
3	Purified Water Purified Water	IBC00016A0	PT Mersifarma
4	Isopropil Alkohol Isopropil Alkohol	A-1087	Smartlab
5	Fluphenazine Hydrochloride RS Fluphenazine Hydrochloride RS	32379-74S	Kanto
6	Fluphenazine Decanoate Dihydrochloride RS Fluphenazine Decanoate Dihydrochloride RS	1280803	USP
7	Fluphenazine Decanoate Injeksi 25 mg/mL Fluphenazine Decanoate Injection 25 mg/mL	IPJ00122	PT Mersifarma

7.2.1.3. Metode yang digunakan : High Performance Liquid Chromatography Method Used

7.2.2. Kondisi pengerjaan HPLC HPLC Operating Conditions

- 7.2.2.1. Kolom : Kolom Brand Ultisil C18 3.9 x 150 x 5 µm (L1) katalog nomor : 5020-01801
Column : Column Brand Ultisil C18 3.9 x 150 x 5 µm (L1) catalog number : 5020-01801
- 7.2.2.2. Aliran pompa : 2,0 mL / menit
Flow rate : 2,0 mL / minute
- 7.2.2.3. Panjang gelombang detektor : 254 nm
Detector wavelengths : 254 nm
- 7.2.2.4. Volume penyuntikkan : 20 µL
Injection volume : 20 µL
- 7.2.2.5. Pelarut : Asetonitril : Metanol : 0,05 M Ammonium Acetate (2:2:1)
Solvent : Acetonitrile : Methanol : 0,05 M Ammonium Acetate (2:2:1)
- 7.2.2.6. Fase gerak : Asetonitril : metanol (50 : 50)
Mobile Phase : Acetonitrile : methanol (50 : 50)
- 7.2.2.7. Waktu Retensi
Retention Time
- a. Fluphenazine Hydrochloride : 2,0 menit ± 1 menit
Fluphenazine Hydrochloride : 2,0 minutes ± 1 minutes
- b. Fluphenazine Decanoate Dihydrochloride : 52 menit ± 10 menit
Fluphenazine Decanoate Dihydrochloride : 52 minutes ± 10 minutes

Good Quality Improves Life's Quality

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Preparation of 0.05 M Ammonium Acetate Solution

- a. Timbang 3,85 gram Ammonium acetate, Masukkan kedalam gelas piala 1000 mL
Weigh 3.85 grams of Ammonium acetate, put it into a 1000 mL beaker.
- b. Tambahkan Purified water hingga 1000 mL aduk dengan magnetic Stirer selama 10 menit hingga larut
Add Purified Water to a volume of 1000 mL, then stir using a magnetic stirrer for 10 minutes until completely dissolved.
- c. Saring larutan menggunakan saringan vakum
Filter the solution using a vacuum filter
- d. Masukkan kedalam botol MPS, sonikasi selama 10 menit.
Put into MPS bottle, sonicate for 10 minutes.

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Preparation of System Suitability Test Solution

7.2.5.1. Larutan Induk UKS Fluphenazine Decanoate Dihydrochloride RS (Konsentrasi : 1 mg/mL

UKS Fluphenazine Decanoate Dihydrochloride RS Main Stock Solution (Concentration: 1 mg/mL)

- Timbang 5 mg Fluphenazine Decanoate Dihydrochloride, masukan kedalam labu ukur 5 mL
Weigh 5 mg of Fluphenazine Decanoate Dihydrochloride, put it into a 5 mL measuring flask.
- Tambahkan 2 mL pelarut, Sonikasi selama 10 menit
Add 2 mL Solvent, Sonicate for 10 minutes
- Tambahkan pelarut sampai tanda batas, kocok sampai homogen.
Add solvent up to the volume, shake until homogeneous.

7.2.5.2. Larutan Induk UKS Fluphenazine Hydrochloride RS (Konsentrasi : 0,8 mg/mL

UKS Fluphenazine Hydrochloride RS Main Stock Solution (Concentration: 0,8 mg/mL)

- Timbang 4 mg Fluphenazine Hydrochloride, masukan kedalam labu ukur 5 mL
Weigh 4 mg of Fluphenazine Hydrochloride, put it into a 5 mL measuring flask.
- Tambahkan 2 mL pelarut, Sonikasi selama 10 menit
Add 2 mL Solvent, Sonicate for 10 minutes
- Tambahkan pelarut sampai tanda batas, kocok sampai homogen.
Add solvent up to the volume, shake until homogeneous.

 Good Quality Improves Life's Quality PT. MERSIFARMA TM. SUKABUMI-INDONESIA	PROTOKOL TRANSFER METODE ANALISIS IDENTIFIKASI, KADAR, BATAS FLUPHENAZINE DAN BATAS CEMARAN YANG TERELUSI LAMBAT DALAM PRODUK JADI FLUPHENAZINE DECANOATE INJEKSI 25 mg/mL TRANSFER PROTOCOL OF ANALYSIS METHOD OF IDENTIFICATION, ASSAY, LIMIT OF FLUPHENAZINE AND LIMIT OF SLOWLY ELUTERATED IMPURITIES IN FINISHED PRODUCT FLUPHENAZINE DECANOATE INJECTION 25 mg/mL		Halaman : 12 dari 15 Pages : 12 of 15
	No. Dokumen	: RD.02-PTP/005.00	
	Document No.	:	
	Tanggal Dokumen Disetujui	:	
	Document Approved Date	:	
	Mengantikan No.	:	
	Supersedes Number	:	

591,87 = Berat molekul Fluphenazine Decanoate

Molecular weight of Fluphenazine Decanoate

664,69 = Berat molekul Fluphenazine Decanoate Dihydrochloride

Molecular weight of Fluphenazine Decanoate Dihydrochloride

7.3. BATAS FLUPHENAZINE DAN BATAS CEMARAN YANG TERELUSI LAMBAT FLUPHENAZINE LIMITS AND SLOWLY ELUTERATING IMPURITIES LIMITS

7.3.1. Persyaratan : Sesuai dengan spesifikasi produk yang berlaku.
Requirement In accordance with applicable product specifications.

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This test does not require re-preparation of the solutions. The analysis is performed using the same HPLC conditions, column, standard solution preparation, system suitability solution preparation, and sample solution preparation as described in the Assay. The percentage of fluphenazine in the Injection is calculated using the chromatograms obtained from the Assay.

7.3.3. Perhitungan kadar Batas Fluphenazine

Limit of Fluphenazine Calculation

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Calculate the recovery (%) limit of Fluphenazine in Fluphenazine Decanoate 25 mg/mL injection solution using the formula:

$$\% \text{ Fluphenazine} = \frac{R_{std} \times W_{std} \times F \times FP_{spi} \times 473,53}{R_{spi} \times V_{spi} \times \text{content} \times FP_{std} \times 510,44}$$

Keterangan :

Description :

R_{spi} = Respons puncak larutan sampel

Peak response of sample solution

R_{std} = Respons puncak larutan standar

Peak response of standard solution

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Fluphenazine Decanoat content per Injection (mg)

FP_{spi} = Faktor pengenceran larutan sampel (250)

Dilution factor of sample solution (250)

FP_{std} = Faktor pengenceran larutan standar (1250)

Dilution factor of standard solution (1250)

F = Kadar Reference Standard (%) sebagai as is

Reference Standard Level (%) as is

473,53 = Berat molekul Fluphenazine

Molecular weight of Fluphenazine

510,44 = Berat molekul Fluphenazine Hydrochloride

Molecular weight of Fluphenazine Hydrochloride