ANALISA PERMINTAAN BARANG

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7.2.3.2. Pembuatan larutan Fase gerak

Preparation of the mobile phase solution

- a. Ukur 920 mL larutan dapar menggunakan gelas ukur, kemudian masukkan ke dalam beaker glass 1 L.
Measure 920 mL of the buffer solution using a graduated cylinder, then transfer it into a 1 L beaker glass
- b. Tambahkan 80 mL asetonitril ke dalam beaker glass, tutup beaker, kemudian aduk menggunakan stirrer selama 15 menit.
Add 80 mL of acetonitrile to the beaker glass, cover the beaker, and stir using a stirrer for 15 minutes.
- c. Pindahkan larutan ke dalam botol MPS, kemudian sonikasi selama 15 menit.
Transfer the solution into the MPS bottle, then sonicate for 15 minutes.

Good Quality Improves Life's Quality
PT. MERSIFARMA TM
SUKABUMI-INDONESIA

PROTOKOL TRANSFER METODE ANALISIS IDENTIFIKASI, KADAR, DISOLUSI, DAN CEMARAN ORGANIK LEVETIRACETAM DALAM PRODUK JADI LEVETIRACETAM 500 MG KAPLET SALUT SELAPUT
ANALYTICAL METHOD TRANSFER PROTOCOL FOR THE IDENTIFICATION, ASSAY, DISSOLUTION TEST AND ORGANIC CONTAMINANTS OF LEVETIRACETAM IN FINISHED PRODUCT LEVETIRACETAM 500 MG FILM-COATED TABLETS

Halaman : 12 dari 22
Pages : 12 of 22
No. Dokumen : RD.02-PTP/103.00
Document No.
Tanggal Dokumen Disetujui :
Document Approved Date :
Menggantikan No. :
Supersedes Number :

7.3. UJI DISOLUSI LEVETIRACETAM SECARA HPLC
DISSOLUTION TEST OF LEVETIRACETAM BY HPLC

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Notes : Perform dissolution tests in accordance with the dissolution equipment usage procedure.

7.3.3. Kondisi pengerjaan HPLC


HPLC Operating Conditions

- 7.3.3.1. Kolom : Inertsil ODS-3V GL-Science C18 150 mm x 4,6 mm 5µm (L1)
Column
Katalog Nomor : 5020-01801
Inertsil ODS-3V GL-Science C18 150 mm x 4,6 mm 5µm (L1) Catalog : Number : 5020-01801
7.3.3.2. Aliran pompa : 1,2 mL / menit
Flow rate
1,2 mL / minute
7.3.3.3. Panjang gelombang detektor : 220 nm
Detector wavelengths
220 nm
7.3.3.4. Volume penyuntikkan : 10 µL
Injection volume
10 µL
7.3.3.5. Fase Gerak : Asetonitril : Dapar pH 5,6 (15 : 85)
Mobile Phase
Acetonitrile : Dapar pH 5,6 (15 : 85)
7.3.3.6. Pelarut : Purified Water
Solvent
Purified Water
7.3.3.7. Waktu Retensi : 2,6 Menit ± 1 Menit
Retention Time
2,6 Minutes ± 1 Minutes
7.3.3.8. Waktu Berhenti : 4 Menit
End Time
4 Minutes

7.3.4. Fase Gerak : Larutan Dapar pH 5,6

Mobile Phase : pH 5,6 Buffer Solution

- 7.3.4.1. Timbang 6,8 gram kalium fosfat monobasa, Masukkan kedalam gelas piala 1000
Weigh 6.8 grams of monobasic potassium phosphate, put it into a 1000 ml beaker.
7.3.4.2. Tambahkan Purified water hingga 1000 mL aduk dengan stirrer selama 15 menit
Add purified water up to 1000 mL and stir with a stirrer for 15 minutes.
7.3.4.3. Atur hingga pH 5,6 ± 0,05 dengan menggunakan KOH encer (0,5 gram KOH dalam 100 mL Purified water)
Adjust to pH 5,6 ± 0,05 using dilute KOH (0.5 grams of KOH in 100 mL of purified water)
7.3.4.4. Saring larutan menggunakan saringan vakum
Filter the solution using a vacuum filter
7.3.4.5. Masukkan kedalam botol MPS, sonikasi selama 10 menit
Put into MPS bottle, sonicate for 10 minutes.

Good Quality Improves Life's Quality

PT. MERSIFARMA TM.
SUKABUMI-INDONESIA

**PROTOKOL TRANSFER METODE ANALISIS
IDENTIFIKASI, KADAR, DISOLUSI, DAN CEMARAN
ORGANIK LEVETIRACETAM DALAM PRODUK JADI
LEVETIRACETAM 500 MG KAPLET SALUT SELAPUT
ANALYTICAL METHOD TRANSFER PROTOCOL FOR
THE IDENTIFICATION, ASSAY, DISSOLUTION TEST
AND ORGANIC CONTAMINANTS OF LEVETIRACETAM
IN FINISHED PRODUCT LEVETIRACETAM 500 MG
FILM-COATED TABLETS**

Halaman	: 15 dari 22
Pages	: 15 of 22
No. Dokumen	: RD.02-PTP/103.00
Document No.	
Tanggal Dokumen Disetujui	:
Document Approved Date	:
Menggantikan No.	:
Supersedes Number	:

7.3.5.4. Perhitungan kadar

Assay Calculation

Hitung kadar (%) Leveteracetam 500 mg dalam kaplet/kaplet salut selaput Leveteracetam 500 mg dengan rumus:

Calculate the assay (%) of Leveteracetam 500 mg in film-coated caplets/caplets using the formula:

$$\% \text{ Leveteracetam} = \frac{R_{\text{spl}} \times W_{\text{std}} \times FP_{\text{spl}}}{R_{\text{std}} \times 500 \times FP_{\text{std}}} \times F$$

Keterangan :

Description :

R_{spl} = Respons puncak larutan sampel
Peak response of sample solution

R_{std} = Respons puncak larutan standar
Peak response of standard solution

W_{std} = Bobot standar yang ditimbang (mg)
Standard weight weighed (mg)

500 = Kandungan Leveteracetam per kaplet atau kaplet salut selaput (mg)
Levetiracetam content per caplets or film-coated caplets (mg)

FP_{spl} = Faktor pengenceran larutan sampel (900)
Dilution factor of sample solution (900)

FP_{std} = Faktor pengenceran larutan standar (100)
Dilution factor of standard solution (100)

F = Kadar Working Standard (%) as is
Working Standard (%) as is

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Phone +62 21 38834021

EXW Xi'an, China

No.	Cat#	Description	Unit Price	Qty	Total
1	13021860	HPLC Column,C18-Plus,4µm,130 A,4.6x250mm,1pc/pk HS CODE: 9027.90	¥3,059.00	2 packs	¥6,118.00
Subtotal					¥6,118.00
Freight Cost to door by UPS					¥298.00
TOTAL					¥6,416.00

NOTE

LEAD TIME: 6-8 Workdays after Receiving the Payment.

DELIVERY: 6-7 Workdays.

PAYMENT: 100% Paid in Advance.

BANK INFORMATION:

Bank Name: CHINA MERCHANTS BANK

Bank Account #: 1299 0556 8210 402

Beneficiary Name: HAWACH SCIENTIFIC CO., LTD

Beneficiary Address: Room 2408, Building A1, No.32 JinYE Rd, Xi'an, Shaanxi Province, P.R.China.

SWIFT Code: CMBCCNBS 231

Bank Branch Address: NO.1 Zhangbayilu, High-tech Development zone, Yanta District Xi'an, Shaanxi Province, China.

PS: Please strictly use the complete Beneficiary Name 'Hawach Scientific Co., Ltd' when making payment.
Any incomplete name may cause failure in receiving your payment.

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Notes :

- Franco Depok
- Ready Stock Before Sold
- Term of payment 100% before delivery

Best Regards,
PT Mass Saintifika Indonesia

Ragita Anggraini
Sales Engineer Executive
0896-0271-9981
ragita@masssaintifikaIndonesia.com