

ANALISA PERMINTAAN BARANG

No	Kode Barang	Nama Barang	Satuan	Stok Saat Ini (A)	Rencana Kebutuhan (B)	Kekurangan (A -B + C)	Pending Penerimaan (C)	Jumlah	Realisasi Pembelian	Realisasi Penggunaan
1	IMO03038	Cotton Swab Sterile	Pc	0	2000	-2000	0	2000	0	0

Catatan :

Swab Viscose Head Polystyrene Stick, 2.2 Mm X 150Mm, 500/Pack, Cat. No : 2122-0009

- Digunakan saat pengambilan sampel mikrobiologi untuk pemeriksaan pada mesin produksi yang sulit dijangkau dengan metode kontak langsung.
- Packing original : 1 Pack @500 Pc.
- Estimasi harga : Rp.592.000/Pack.


Latar Belakang Permintaan

- Dalam rangka mendukung pelaksanaan cleaning validation diperlukan pengambilan sampel mikrobiologi dengan cara swab untuk memastikan bahwa proses pembersihan mesin produksi mampu menghilangkan mikroorganisme Sehingga diperlukan pengadaan cotton swab steril sebagai bahan pendukung kegiatan pengambilan sampel.
- Masa frekuensi pakai Cotton Swab Sterile yaitu setiap kali validasi pembersihan.
- Pemesanan pertama kali Dept. QC. Terakhir pemesanan melalui PPIC pada bulan Agustus 2024 (PR151675)

Terlampir rencana kebutuhan dan penawaran harga.

Dibuat

Mengetahui

 Dra. Suci
15/09/26

 Dra. Ayu
15/09/26



PT INTRALAB EKATAMA

Jl. Terapi Raya AD 2 - Bumi Menteng Asri
Bogor, Indonesia, 16111
Phone : 0251-8359110, 8311662
E-mail : intralab@indo.net.id, info@intralab.co.id



QUOTATION

To :

PT MERSIFARMA TIRMAKU MERCUSANA

Attn : Bpk. Naufal Hisyam Fasya

Quo No. : FYAP/0215/VIII/24

Date : 16 Agustus 2024

Payment term : net 30

No.	Catalog No.	Description	Quantity	Unit Price	Disc	Total	Tax
1.	2122-0009	Swab Viscose Head Polystyrene Stick, 2.2Mm X 150Mm , 500/Pack Note : Kemasan 500 pcs/pack, READY STOCK	1 PACK	Rp 592,000	%	Rp 592,000	11%
2.	2122-0022	Swab Cotton Head Polystyrene Stick 2.2Mm X 150Mm , 4/Case Note : Kemasan 500 pcs/pack, INDENT 4-6 BLN	1 PACK	Rp 592,000	%	Rp 592,000	11%

*) Untuk barang ready, stock tidak mengikat / stock berjalan.

Notes :

Sub Total	: Rp	1.184.000
PPN (11%)	: Rp	130.240
Other	: Rp	0
Total	: Rp	1.314.240

Bank Account :

PT BANK SYARIAH INDONESIA(BSI)

KCP BOGOR TAJUR 1

A/N. : PT INTRALAB EKATAMA

A/C. : 105.452.2165

BANK CODE : 451

Best Regards,

Frans Yogas Aria Putra

*Jika
pasihan CoA terdapat jika tidak
ganti supplier yg tersedia*

No	Line	Nama Mesin	Produk Marker	Juli	Agust	Sept	Okt	Nov	Des
No.	Line	Machine's Name	Marker Product	July	August	Sept	Oct	Nov	Dec
Unit 1 Non Steril									
1	Suspensi	Peralatan Penimbangan	Corimoxazole 400/40 mg suspensi				4	4	4
		High Shear Emulsifying Machine					6	6	6
		Mixing Tank 2000 L (double jacket)					6	6	6
		Filling & Capping Shanghai Xiongshun					12	12	12
2	Sirup	Peralatan Penimbangan	Paracetamol 60 mL Syrup				6	6	6
		Turbo Mixer					3	3	3
		Mixing Tank 2000 L					5	5	5
		Holding Tank 2000 L					6	6	6
		Filling & Capping Shanghai Xiongshun					12	12	12
3	Tablet (uncoated)	Peralatan Penimbangan	Cotrimoxazole 400/80 Tablet			6	6	6	
		OG-YK 160				7	7	7	
		FBG FL-250				10	10	10	
		Automatic Lifting Bin Mixer				3	3	3	
		Accura Press IV 61 STN				37	37	37	
		Blister Accurate 240 Plus				7	7	7	
Unit 2									
5	Kaplet	Peralatan Penimbangan	Amoxicillin Trihydrate 500 Kaplet		4	4	4		
		Vacuum Transfer			3	3	3		
		V- Mixer 2500 L			2	2	2		
		Fluidpack Accura Press IV 45			32	32	32		
		Stripping Chunguang DXDP 350			7	7	7		
Unit 3									
6	Tablet (uncoated)	Peralatan Penimbangan	Allopurinol 100 Tab			4	4	4	
		OG-YK 160				7	7	7	
		Super Mixer GHL 1200				6	6	6	
		FBG FL-250				6	6	6	
		Bin Blender HZD 2000				10	10	10	
		Turbo Mixer BRT				3	3	3	
		Fluidpack Accura Press IV 75				42	42	42	
		Blister Accurate 240 Plus				8	8	8	
7	Tablet (coated)	Peralatan Penimbangan	Cetirizine				4	4	4
		OG-YK 160					6	6	6
		Super Mixer GHL 1200					6	6	6
		FBG FL-250					6	6	6
		Bin Blender HZD 2000					10	10	10
		Turbo Mixer BRT					3	3	3
		Fluidpack Accura B4 35					32	32	32
		Coating Changzou BGB 400					8	8	8
		Stripping Chunguang DXDP 350					8	8	8
Unit 4									
9	Tablet (uncoated 2)	Peralatan Penimbangan	Alprazolam 0,5		4	4	4		
		Ball Mill			10	10	10		
		Super Mixer SHK 200			6	6	6		
		FBD FL-60			6	6	6		
		OG-YK 160			6	6	6		
		Fluidpack ATX-1			32	32	32		
		Fluidpack Accura B4 35			32	32	32		
		Blister Accurate 240 Plus			8	8	8		
11	Tablet (coated 2)	Peralatan Penimbangan	Haloperidol 5 TSS				4	4	4
		Ball Mill					10	10	10
		Super Mixer GHL-400					6	6	6
		Planetary Mixer					3	3	3
		FBD FL-120					6	6	6
		OG-YK 160					6	6	6
		Fluidpack Accura 35, 75					42	42	42
		Coating Chengdong, Changzou BGB 400					8	8	8
		Stripping Chunguang DXDP 350					8	8	8

**mersi**PT. MERSIFARMA TM
SUKABUMI-INDONESIA**RENCANA KEBUTUHAN COTTON SWAB STERIL VALIDASI
PEMBERSIHAN 2026**

No	Line	Nama Mesin	Produk Marker	Juli	Agust	Sept	Okt	Nov	Des
No.	Line	Machine's Name	Marker Product	July	August	Sept	Oct	Nov	Dec
12	Kapsul	Peralatan Penimbangan	Fluoxetine HCL 20				4	4	4
		Conical Mixer					4	4	4
		Filling Capsule NJP					14	14	14
		Stripping Chunguang DXDP 350					8	8	8
Unit 1 steril									
14	Ampul (1)	Peralatan Penimbangan	Haloperidol injeksi				4	4	4
		Mixing Tank 50 L					7	7	7
		Transfer Tank 50 L					7	7	7
		Compact Line Filling Ampoule					17	17	17
16	Vial	Peralatan Penimbangan	Levofloxacin Infus				4	4	4
		Mixing Tank 500 L					7	7	7
		Transfer Tank 500 L					7	7	7
		Filling Vial					17	17	17
Unit 6A									
18	Vial	Peralatan Penimbangan	Levofloxacin Infus			4	4	4	
		Mixing Tank 1000 L				7	7	7	
		Transfer Tank 500 L				7	7	7	
		Filling Vial				17	17	17	
Jumlah				0	152	343	679	527	336
Total				2037					

Keterangan :

Catatan Jumlah cotton swab kemungkinan akan berubah
 15/70
109

Citotest Labware Manufacturing Co.,Ltd

Certificate of Analysis

Lot No.:231881

Product Name:Collection Swabs			Part No.:2122-0009				
Raw Material:VISCOSE、PS			Specification :PS stick+VISCOSE head				
Standard:Product inspection criterion&PO							
Package :1pcs/peel-pk, 500pcs/disp.box, 2,000pcs/carton							
Quantity:10,000pcs			Report Date:2024.05				
Test item	Technical requirement	Inspection level	Sample	Qualified quality level(AQL)		Defect Q'ty (pcs)	Inspection result
Dimension	1,Height of the swab : 150±3mm 2,Diameter of the head: 5±0.5mm 3,Length of the viscose head: 18±3mm	5pcs/lot	5pcs	0		0	PASS
				Ac	Re		
				0	1		
Appearance	1,Swab head should be soft ,flexible. 2,No dirty,no odour,no stain in the head.	GB/T2828.1 - 2012/ISO2859-1:1999 Normal test level II	200pcs	1.5		0	PASS
				Ac	Re		
				7	8		
Function	1,The swab head should not drop when exerting the force of 4.9N to pull for 10s. 2,The swab should not drop from the tube when exerting the force of 9.8N to pull for 10s. 3,Soakage not less than100ul/pieces. 4,PH value:the liquor shouldn't show pink when add phenolphthalein in,and shouldn't show yellow if add bromcresol purple in. 5,No blue light under the UV-light. 6.Sterilization effect meet the standard of ISO11137	10pcs/lot	10pcs	0		0	PASS
				Ac	Re		
				0	1		
Package	Quantity:Conform to requirements of PO	10 inner boxes/lot 1 outer box/lot	10 inner boxes 1 outer box	0		0	PASS
	Packing type:Conform to requirements of PO			Ac	Re		
	Label and marking:Clear, correct and of correct location			0	1		
Conclusion	The products meet all the requirement by our test						

Inspected by:Wangjing

Checked by:Wang Jing

COMMERCIAL CENTER&R-D Center

10F, Building 3, 89 Shengji Rd, Jiangning
District, Nanjing 211106, China
Tel.: 0086 25 86216803
Fax: 0086 25 86214463
E-mail: info@citotest.com

PRODUCTION

No. 339 Beihai West Road, Haimen, 226100 Jiangsu, P.R. China
Tel.: 0086 513 82259348
Fax: 0086 513 82110195
E-mail: info@citotest.com.cn

**CITOTEST LABWARE MANUFACTURING CO.,LTD****www.citotest.com****EXECUTIVE SUMMARY**

For MD Registration – Ministry of Health, Republic of Indonesia

Section	Details
1. Device Name	Trade Name: Citolab Model / Type: Citolab Classical Swab Microbiological Collection Swab viscose head + PS stick (2122-0009) Manufacturer: CITOTEST Labware Manufacturing Co.,Ltd Country of Origin: China
2. Device Description	Microbiological collection swab, Dim. 2.2mm(dia)x 150mm (6"), Viscose head with Dia. 5mm approx., PS stick. individual peel pack, GAMMA sterile
3. Marketing History link FSC	First marketed in 2008 in Spain. Now sold in 25+ countries including Singapore, Indonesia, Brazil, and Australia. Certified CE (Class Is, Rule 5) <i>Below is CFSC link</i> https://www.nmpa.gov.cn/datasearch/search-info.html?nmpa=aWQ9ODU5YmQ5MDYwY2FhZGUyOGYwMmY0MDhhYTMwODMzYmMmaXRlbUlK PWZmODA4MDgxODNjYWQ3NTAwMTgzY2I2Nm ZlNjkwMjg1
4. Intended Use	The dry swabs are designed for biological sampling from surface and wide orifices, Such as skin, mouth, nose and throat, for Processing samples and isolation in culture medium, and soon. The swabs series provided with a complete range of sticks and tips including cotton swabs, viscose swabs and polyester swabs.
5. Indications for Use	Peel the dry swabs units; Just hold the swab to sampling from surface and wide orifices; Record the required information and send to the laboratory with minimum delay; The used swabs must be disposed of According to laboratory regulations for infectious waste.
6. Safety Information	Store the swabs at dry and dark place in Room temperature (2 - 30°C). Use before the expiration date which is clearly printed on the product unit.
7. Performance Information	The shelf life of Classical Swab is 3 years. The unopened or undamaged units are guaranteed sterile. Do not use, if there is a breakage of the pouch, a leakage of the tube

COMMERCIAL CENTER&R-D Center

10F, Building 3, 89 Shengli Rd, Jiangning
District, Nanjing 211105, China
Tel.: 0086 25 86216803
Fax: 0086 25 86214463
E-mail: info@citotest.com

PRODUCTION

No. 339 Beihai West Road, Haimen, 226100 Jiangsu, P.R. China
Tel.: 0086 513 82259348
Fax: 0086 513 82110195
E-mail: info@citotest.com.cn



CITOTEST LABWARE MANUFACTURING CO.,LTD

www.citotest.com

Prepared by:

Name : Janna Ren

Position : International Sales Department

Date : 18/12/2025

