



ANALISA RENCANA PEMAKAIAN REAGEN

Tanggal SBP :

No SPB :

1. Anaerogas pack 3,5 L (Thermo Scientific No. Cat 4323LH); (ex Himedia No. Cat LEO02A)

Serapan Material	Jun 26	Jul 26	Aug 26	Sep 26	Oct 26	Nov 26	Des 26
Stok awal	3 pcs						
Aktual pemakaian							
Rencana kebutuhan	5 pack @ 10 pcs (50 pcs)						
Outstanding PO							
Rencana pemakaian							
- Kedatangan Bahan			1 pcs	1 pcs	1 pcs	1 pcs	1 pcs
- Stabilitas Produk Baru		1 pcs	2 pcs	3 pcs	2 pcs	1 pcs	
- Validasi Proses	6 pcs	4 pcs	8 pcs	4 pcs	4 pcs	4 pcs	
Stok sisa	47 pcs	42 pcs	31 pcs	23 pcs	16 pcs	10 pcs	9 pcs

Note : sisa untuk stok reagen tahun 2026

Klaten, 10 Juni 2026

Disusun oleh,

Nilam Wardah

Diperiksa Oleh,

Yeni Tri Hastuti

Paulijanto

Disetujui oleh,

Suherman



Perhitungan :

Juni	Juli	Agustus	September
1. Validasi Proses Gamangin 10 Sampel = 2 pcs 2. Validasi proses OBK 12 sampel = 4 pcs	1. Stabilitas 6 sampel = 1 pcs 2. Validasi proses Phymogala 12 sampel = 4 pcs	1. Kedatangan bahan 5 sampel = 1 pcs 2. Stabilitas 10 sampel = 2 pcs 3. Validasi proses OBK 12 sampel = 4 pcs 4. Validasi proses phymogala 12 sampel = 4 pcs	1. Kedatangan bahan 5 sampel = 1 pcs 2. Stabilitas 16 sampel = 3 pcs 3. Validasi proses gamangin 12 sampel = 4 pcs
Oktober	November	Desember	
1. Kedatangan bahan 5 sampel = 1 pcs 2. Stabilitas 10 sampel = 2 pcs 3. Validasi proses phymofit 12 sampel = 4 pcs	1. Kedatangan bahan 5 sampel = 1 pcs 2. Stabilitas 4 sampel = 1 pcs 3. Validasi proses phymogala 12 sampel = 4 pcs	1. Kedatangan bahan 5 sampel = 1 pcs	